



EmMa.Test

Feline

MAJOR XM + MINOR XM + BLOOD TYPING A+B
RESULT FORM

Date :

Test performed by :

EmMa.Test lot number :

To interpret your result, stick the membrane below, **respecting the arrow color and the orientation of the membrane. Read the result immediately after testing.**

MAJOR XM

Donor ID
(pink membrane) :

Donor
Blood Typing :

A B

AB

Major XM

+ -



MINOR XM

Recipient ID
(grey membrane) :

Recipient
Blood Typing


A B

AB

Minor XM

+ -

COMPATIBLE = SAFE TRANSFUSION
(ONLY 3 POSSIBLE RESULTS)



No line or white line = Negative result

A weak line must be considered as a positive line.