



EmMa.Test

Feline

MAJOR XM + MINOR XM + BLOOD TYPING A+B
RESULT FORM

Date : _____

Test performed by : _____

EmMa.Test lot number : _____

To interpret your result, stick the membrane below, **respecting the arrow color and the orientation of the membrane.** Read the result immediately after testing.

MAJOR XM


Donor ID
(pink membrane) :

Donor
Blood Typing :

A B
AB

Major XM

+ -



A
B
XM



A
B
XM

MINOR XM

Recipient ID
(grey membrane) :


Recipient
Blood Typing

A B
AB

Minor XM

+ -

COMPATIBLE = SAFE TRANSFUSION (ONLY 3 POSSIBLE RESULTS)



Red line = Positive result
No line or white line = Negative result

A weak line must be considered as a positive line.

